

carilex®

# Entrix NX

It is the best example of affordable quality. Entrix NX benefits from Carilex's deep knowledge of the therapeutic support surface. It is equipped with the most sophisticated Carilex IPS technology, which guarantees the best performance for low and medium-risk patients.

The pump is suitable for two different mattress:

- Entrix Mattress Overlay
- Lux Mattress Replacement

Entrix Mattress Overlay

LUX Mattress Replacement



## 1 Intelligent Pressure Sensing

Responds to patient movements on the mattress by automatically adjusting internal mattress pressure.

## 2 Low-Pressure Alarm

Audible and visible alert to notify when abnormal internal pressure occurs.

## 3 Mattress Size Selection

Compatible with multiple mattress sizes, which can be selected at start-up.

## 4 Dynamic / Static Therapy Modes

Provide both dynamic and static functions to meet individual patients' needs.

## 5 Power Failure Alarm

Audible and visual power failure alarms alert caregivers in case of an abnormal power outage.

## 6 Memory Recall

Restores to the previous setting automatically after an abnormal power outage.

## 7 Alarm Mute

Disable the audible indicator system.





## Key Features



Alternating Pressure Relief



Lightweight Power Unit



Pressure Redistribution



Whisper Quiet



Multiple Weight Setting



Memory Recall



Mattress Overlay



Durable Material



Multiple Mattress Sizes



Mattress Replacement



Easy Cleaning



5" Air Cell



Cable Management



Protective PE layer



7" Air Cell

### Technical Data

Power Unit Model	REF : SR366
Power Input	100-120Vac 60Hz 0.2A / 220-240Vac 50/60Hz 0.2A
Power Consumption	Max 15W
Weight	<ul style="list-style-type: none"> <li>Power Unit : 2.8kg (6.2lb)</li> <li>Mattress Set : 6kg, 7kg (13lb, 16lb)</li> </ul>
Dimension	Power Unit : 12" x 9" x 5" (29.5 x 25.5 x 12 cm) ( W x H x D ) Overlay Mattress Set ( W x L x H ) : 32" x 80" x 5" ( 82 x 200 x 13 cm ) 36" x 80" x 5" ( 90 x 200 x 13 cm ) Replacement Mattress Set ( W x L x H ) : 32" x 80" x 7" ( 82 x 200 x 18 cm ) 36" x 80" x 7" ( 90 x 200 x 18 cm ) other size upon request
Agency Approval	<ul style="list-style-type: none"> <li>IEC/EN60601-1</li> <li>ANSI/AAMI ES60601-1</li> <li>AS/NZS IEC60601.1</li> </ul>
Operating Condition	<ul style="list-style-type: none"> <li>Temperature : 5°C to 40°C (41°F to 104°F)</li> <li>Humidity : 15% to 90%</li> </ul>
Transportation & Storage Condition	<ul style="list-style-type: none"> <li>Temperature : -25°C to 70°C (-13°F to 158°F)</li> <li>Humidity : 0% to 90%</li> </ul>
Weight Capacity	<ul style="list-style-type: none"> <li>Entrix : 150kg (330lb)</li> <li>Lux : 180kg (400lb)</li> </ul>

This Power Unit is not AP/APG protected

[www.carilexmedical.com](http://www.carilexmedical.com)

Carilex Medical, Inc. **HQ**

TEL : +886-3-328-7882  
 Email : [amea@carilexmedical.com](mailto:amea@carilexmedical.com)  
 No.77, Keji 1st Rd., Guishan District,  
 Taoyuan City (33383), Taiwan

Carilex Medical USA

TEL : +1-866-551-4227  
 Email : [americas@carilexmedicalusa.com](mailto:americas@carilexmedicalusa.com)  
 17068 Evergreen Place, Unit A, City of Industry,  
 CA 91745 USA

Carilex Medical DE

TEL : +49-4121-476896  
 Email : [eu@carilexmedical.com](mailto:eu@carilexmedical.com)  
 Fritz-Strassmann Str. 5, 25337 Elmshorn  
 Germany

Carilex Medical NL

TEL : +31-2089-05788  
 Email : [nl@carilexmedical.com](mailto:nl@carilexmedical.com)  
 Overschiestraat 59, 1062XD  
 Amsterdam

Carilex Medical UK

TEL : +44-1323-761900  
 Email : [uk@carilexmedical.com](mailto:uk@carilexmedical.com)  
 5th Floor, 167-169 Great Portland St,  
 London W1W 5PF, United Kingdom