



## Kangaroo™ 924 Enteral Feeding Pump

### FEATURES

- Accurate pre-set volumetric delivery with volume totalizer
- Intuitive and user-friendly controls
- Easy-to-read LED display
- Built-in pole clamp
- Over-infusion safeguard
- Audio and visual alarms
- 24-hour memory
- DEHP-free pump sets

### THE STANDARD IN ENTERAL FEEDING

When ease of use, accuracy and dependability matter, the Kangaroo 924 Enteral Feeding Pump delivers it all.

### EASY-TO-USE CONTROLS

- Makes set-up and operation intuitive.

### PRESET VOLUME AND OVER-INFUSION SAFEGUARD

- Ensures correct and consistent delivery of nutrition.

### AUDIBLE AND VISUAL ALARMS

- Reinforce clinician confidence and patient safety.

### HIGHLY VISIBLE LED DISPLAY AND BUILT-IN POLE CLAMP

- Aids in ease of use.



Easy-to-use pump set



1000 ml pump set



924 spike set



Kangaroo safety screw spike

## FEATURES

- Pumping Mechanism – Rotary peristaltic
- Operating Range – 1 - 300 ml/hr in 1 ml/hr increments
- Volume to be Delivered – 1-2000 ml in 1 ml increments
- Accuracy – +/- 10%
- Occlusion Pressure – 12 psi
- Battery – Nickel Metal Hydride, Rechargeable
- Battery Life – 16 hours at 125 ml/hr
- Memory – 24 hours
- Dimensions – 7.25" x 6.6" x 4.33"
- Weight – 3.89 lbs.
- Alarms – 6 options
- Multilingual – 13 languages

## ORDERING INFORMATION

Code Number	Description	Case Quantity
392455	Kangaroo 924 Enteral Feeding Pump with pole clamp	1
773621	1000 ml Bag Set (Indiv. Pack)	30
703621	1000 ml Bag Set (Bulk Pack)	30
774625	Spike Set (Indiv. Pack)	30
772025	500 ml Bag Set	30

### KANGAROO FEEDING SETS FOR USE WITH NESTLÉ™\* SPIKERIGHT™\* SYSTEM

775756	Kangaroo Safety Screw Spike for Kangaroo 924 Pump	30
775762	Safety Screw Spike Adapter Cap	30

[www.kangarooepump.com](http://www.kangarooepump.com)



™\* Trademark of the Société des Produits Nestlé S.A.

COVIDIEN, COVIDIEN with logo and ™ marked brands are trademarks of Covidien AG or an affiliate. © 2008 Covidien AG or an affiliate. All rights reserved.

H6266-5M-0408